

You can find out everything about the successful SQOOM concept at [www.sqoom.com](http://www.sqoom.com)

### SQOOM clinical


110.201201.EN

A registered trademark of Schick Medical GmbH

Auf dem Wall 29  
D-78628 Rottweil

Phone +49(0)741/1741240  
Fax +49(0)741/17412419  
info@SQOOM.com  
www.SQOOM.com

SQOOM clinical is for use in medical practices and clinics. Request our detailed documents. You can receive further information about the SQOOM concept from your doctor or pharmacist.



### The new clinically tested anti-aging concept.

SQOOM is an innovative anti-aging treatment that regenerates your skin from deep within. While the SQOOM device is responsible for the transport, the SQOOM gels ensure the provision of the right and active agents.



### The innovative SQOOM concept for skin problems and sensitive skin.



the concept of beauty

the medConcept  
[www.SQOOM.com](http://www.SQOOM.com)

## SQOOM acts deep inside the skin

In contrast to conventional creams and gels, which are rubbed into the skin by fingertips or with the palm of the hand, the SQOOM concept channels highly effective cosmetics deep into the skin - reaching tissue layers that the active agents of conventional creams cannot.


SQOOM transports valuable active agents into the deeper skin and tissue layers, using a transducer that generates one million soft vibrations per second. Innovative biomolecule preparations from thermal, mechanical and physiochemical components combined with the microcomputer-controlled SQOOM enable a highly complex effect.

SQOOM is an innovative microcomputer-controlled treatment system using highly developed active gels to regenerate the skin from deep within. The worldwide patented „Synchron-Technology“ features a unique combination of iontophoresis with ultrasound and active cosmetics. The gentle treatment process regenerates the connective tissue, effectively improving the skin structure and relieves pain.



## The valuable SQOOM gels

The worldwide patented IONOZYM system features a unique combination of iontophoresis and ultrasound (sonophoresis). In other words: the ultrasound is the door bell while the ionisation subsequently opens countless little doors through which the highly effective agents can penetrate into the skin - millions of times per second.



Impure, oily skin as typically seen in cases of acne.



Scars, as typically related to surgery and injuries.



For a targeted massage of back strain and muscle tension.



Dry, irritated skin, as typically related to neurodermatitis.

Dry, flaked skin, as typically related to psoriasis.



Reddened skin, as typically related to rosacea and facial erythrosis.

